
synapsis® point 2

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ENGLISH

Table of contents

1.	General Description, Purpose and Effect	22
1.1	General Description	22
1.2	Medical Purpose	22
	Literature about the electro-therapy stimulation (induction therapy)	22
1.3	User Groups, Prospective patients, Environment	22
2.	First Use - Proper Handling	23
2.2	This is how you start a treatment	23
3.	Contraindications, Precautions and Adverse Effects	24
3.1	Contraindications	24
3.2	Side-effects	24
3.3	Precautions	24
4.	Cleaning and Maintenance	25
5.	Inspection, Maintenance and Repairs	25
5.1	State of the Battery Charge	25
5.2	Monitoring Proper Function	25
5.3	Maintenance	26
6.	Warranty	26
7.	Technical Data of the Control Unit	26
8.	Applied Standards and Guidelines	27
	Appendix/Symbols:	28
	Appendix 1 – EMV-Tables	31

1. General Description, Purpose and Effect

1.1 General Description

The electro-therapy system Synapsis Point 2 produces low frequency pulsations of very low voltage amplitude and specific frequency scanning. They are administered at the wrists or at selected cutaneous zones (stratum corneum) and acupuncture points.

As very low voltages are used, only an infinitesimal energy transfer takes place. The effect results from the resonance of the applied rhythmic signals with the sensitive sensory mechanisms of the skin surface and their regulatory interaction with functionally important cerebral areas. The administered specific frequency patterns encompass the range from 1 Hz to 100 Hz and are thus corresponding to the frequency sequences of the human brain.

1.2 Medical Purpose

The purpose of the generic Synapsis family of products is to support the treatment of depression, mental exhaustion/stress and insomnia through electro-therapy stimulation. The application is through cuffs at the wrists and/or sticky electrodes or the use of a point applicator at zones and acupuncture points of the skin (stratum corneum).

The International Mandel Institute, Esogetics GmbH and partners of the Esogetics GmbH are offering trainings and courses for the electro-therapy (induction therapy).

Literature about the electro-therapy stimulation (induction therapy):

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1.3 User Groups, Prospective patients, Environment

The intended users are instructed trained therapists or professionals (therapist, physician, naturopath), as well as patients who have been instructed by a therapist.

The Synapsis Point 2 instrument is used in the clinical setting.

The Synapsis Point 2 instruments can be used with all patients starting at the age of 9, independent of size, weight or gender.

2. First Use - Proper Handling

illustr. 1 – Instrument, point applicator, manual and case

The electro-therapy system Synapsis Point 2 consists of a convenient battery-powered regulatory unit, a point applicator and a handheld electrode, as well as the manual. The instrument and its accessories are kept in a plastic case (illustr. 1).

illustr. 2 – Instrument/ socket

The Synapsis Point 2 instrument has a socket for the hookup of the connector cable (illustr. 2).

The connector cable of the Synapsis Point 2 instrument possesses an interface for the attachment of the point applicator accessory.

illustr. 3 – Instrument and point applicator

To get the instrument ready, please insert the round plug of the point applicator into the socket on the top of the instrument, aligning the two red dots (illustr. 2 and 3). The instrument is now ready for use.


2.2 This is how you start a treatment

illustr. 4 – Synapsis Point 2 device with buttons from above

By pushing the ① button the device is switched on. While the patient is holding the cylindrical electrode in his hand, the practitioner is positioning the point applicator on a selected zone or acupuncture point (illustr. 5) of the skin (stratum corneum). Then the program button Beta, Delta, Alpha, Theta or Gamma is pushed to start the chosen program (illustr. 4).

The running application is visualized through a blinking indicator light. A repeated pushing of the program button Beta, Delta, Alpha, Theta or Gamma will terminate the program immediately at any

illustr. 5 – Example of use

time. After the completion of the program sequence the instrument switches off automatically. The round plug can remain in the socket in between applications. If you are not using the device for an extended period, you might want to remove the round plug.

3. Contraindications, Precautions and Adverse Effects

3.1 Contraindications

- Pregnancy
- Severe myocardial disease or severe cardiac arrhythmia
- Cardiac pacemaker, especially with treatments in the upper part of the body
- Treatment near the carotid sinus (on the side of the neck)
- Treatments on damaged skin, or open wounds, or other invasive treatment.

3.2 Side-effects

In adherence to the above listed contraindications and the designated use of the instrument there are no known side-effects. In case there are any undesirable side-effects or events, they should immediately be reported to the manufacturer Esogetics GmbH in Bruchsal.

3.3 Precautions

In case of any of the following symptoms or indications it is mandatory to consult your physician or non-medical practitioner prior to treatment with the Synapsis Point 2:

- highly inflammatory illnesses accompanied by fever,
- advanced, malignant tumors,
- thromboses,
- psychoses,
- epilepsy,
- spastic paralysis,
- inflammatory skin conditions,
- metallic implants in the area of the treatment,
- active implants.

The Synapsis Point 2 device may never be used in the vicinity of a therapeutic shortwave or microwave device (minimum distance requirement 1 m = 3 ft.) The simultaneous treatment with the Synapsis Point 2 System and a high-frequency surgical instrument is inadmissible.

4. Cleaning and Maintenance

No special maintenance is required for the Synopsis Point 2. You can clean the casing and accessories with a damp cloth and regular, mild household cleansers and disinfectants.

Attention: Never clean the device with organic solvents or abrasive cleansers. These substances may damage the casing or the keypad.

5. Inspection, Maintenance and Repairs

5.1 State of the Battery Charge

If battery power is too low, a “Low-Battery” message will appear on the display. To change the battery, open the battery case at the back of the instrument and remove the battery

Attention: Batteries are hazardous waste and must be disposed of according to regulations. (Specialty retailers usually supply battery disposal boxes.) Observing the correct polarity, attach the battery clip to the new battery and close the battery case.

Caution: Disconnect the patient from the device prior to battery replacement. If you are not using your equipment for an extended time, please remove the battery. Any battery leakage might otherwise destroy your unit.

5.2 Monitoring Proper Function

illustr. 6 – Device/ socket/ indicator bulb

In order to monitor its proper function the instrument has an indicator bulb on top of the instrument beside the port for the round plug (illustr. 6).

To test the proper function, please insert the round plug into the port on top (illustr. 6). Attach the point applicator and start an application, as described in chapter 2.2

illustr. 7 – Point applicator and hand electrode

Touch the point applicator (illustr. 7) to the hand electrode (illustr. 7). In case of proper function the orange indicator bulb (illustr. 6) will be flashing. Terminate the program by pushing the program button once more.

5.3 Maintenance

There is no special maintenance required for the instrument. In case of uncertainties or complaints about the functionality, please contact the manufacturer directly: Esogetics GmbH in Bruchsal, Germany.

Only the manufacturer or persons appointed by the manufacturer are authorized to repair or make changes to the instrument.

We assume no liability for inappropriate changes and handling that does not conform to the described purpose of use.

6. Warranty

The warranty period for the instrument is one year from the date of purchase. Please contact the manufacturer directly, if your instrument does not work properly although you have handled it correctly. Your faulty unit will be immediately repaired or substituted.

7. Technical Data of the Control Unit

Type designation:	Synapsis Point 2
Power supply:	9 V E-Block / 6LR61 / 6LF22
Nominal current:	17 mA
Type of construction:	Hand-held unit
Protection rating:	III
Humidity protection rating:	IP 65
Duration of treatment:	30 second timer with auto-off 30 sec. after the last keystroke
Auto-Off:	At the end of the treatment
Output voltage (unloaded 1 kOhm)	5V pp; 1,6 V rms
Output voltage (loaded 1 kOhm)	5V pp; 0,5 V rms
Output current (loaded 1 kOhm):	0,5 mA
Output frequency:	1 Hz – 100 Hz
Output current (loaded 1 kOhm)	0,5 mA
Weight:	155 g (including battery)
Measurements:	141 x 63 x 34mm
Service conditions:	+5°C to 35°C
Storage:	-20°C to 70°C
Air pressure:	800 hPa to 1030 hPa
Air humidity:	30% to 93%

Contents of the package:

Control unit Synopsis Point 2, instruction manual, point applicator and a hand electrode.

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Classification according to the Medical Device Directive Class, Appendix 9, Rule 9, Paragraph 1: Class IIa:

**Manufacturer:**

Esogetics GmbH, Hildastr. 8, D-76646 Bruchsal, Germany

Telefon: +49 (0) 7251/800 10, Fax: +49 (0) 7251/800 155

Web: www.esogetics.com

Mail: info-de@esogetics.com

Download of the instruction manuals under:

www.esogetics.com/downloads

Esogetics is a registered medical device manufacturer. We adhere to a certified and monitored quality management system according to Medical Device Directive 93/42/EWG and to DIN EN ISO 13485 with an implemented monitoring and reporting system.

Appendix/Symbols:

The conformity of the device with the relevant EU directive 93/42/EWG MDD is validated through the CE sign with the number of the identified post.

Refers to the necessity for the user to consult the manual

Device type BF according to IEC/EN 60601-1

Refers to the necessity for the therapist to consult the manual in regard to important security related statements, like warnings or cautionary alerts, which could not be attached to the medical product directly for a variety of reasons.

Indicates the manufacturer of the medical product according to the EU guidelines 90/385/EWG, 93/42/EWG and 98/79.

The crossed out trashcan symbol on your product signifies that this product is an electrical and electronic device that is subject to special waste disposal regulations. To enforce recycling, and WEEE waste management (Waste Electrical and Electronic Equipment) and to save the environment and protect health, the European regulations proposes two options for the selective collection of devices taken out of service:

- The dealer takes the old item back, when you buy a new one.
- Old items can be recycled at specific collection sites.

Defines the temperature limits that the medical product can safely be subjected to.

Defines the range of humidity in which the medical product functions securely.

Defines the air pressure, which the medical product safely tolerates.

Serial Number

Serial Number

The Serial no. of your Synapsis device is configured as follows:

XX | **JJJJ** - **xxxxxx** — Continuous serial no.

└─ Production year
└─ Device Type

(**SP 2** = „point 2“; **SW 2** = „wave 2“; **SH 2** = home 2“)

Example:

SH 2 2017 03937 = synopsis home 2 from 2017 with the no. 3937.

Appendix 1

EMV Tables

**Guidance and manufacturer's declaration -
electromagnetic emissions**

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the synopsis point 2, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The synopsis point 2 is intended for use in the electromagnetic environment specified below. The customer or the user of the synopsis point 2 should assure that it is used in such an environment.

Not mentioned test sections are not applicable for synopsis point 2.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The synopsis point 2 uses RF energy only for its internal function. Therefore, RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The synopsis point 2 is suitable in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration - electromagnetic immunity

The synopsis point 2 is intended for use in the electromagnetic environment specified below. The customer or the user of the synopsis point 2 should assure that it is used in such an environment (point healthcare environment – excluding outdoors, vehicles and public places).

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV ± 15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
magnetic field Power frequency (50/60) Hz magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	$V_1 = 6$ Vrms 150 kHz – 80 MHz $V_1 =$	6 Vrms	
Radiated RF IEC 61000-4-3	$E_1 = 10$ V/m for point healthcare environment 80 MHz – 2,7 GHz	10 V/m	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment

Guidance and manufacturer's declaration - high frequency wireless communication devices

Test frequency (MHz)	Service	Maximum power (W)	Distance (m)	Test level (V/m)
385	TETRA 400 (380...390MHz)	1,8	0,3	27
450	GMRS 460, FRS 460 (430...470MHz)	5	0,3	28
710, 745, 780	LTE Band 13/ 17 (704...787MHz)	0,2	0,3	9
810, 870, 930	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5 (800...960MHz)	2	0,3	28
1720, 1845, 1970	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1/ 3/ 4/ 25, UMTS (1.7...1.99GHz)	2	0,3	28
2450	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band7 (2.4...2.57GHz)	2	0,3	28
5240, 5500, 5785	WLAN 802.11 a/ n (5.1...5.8GHz)	0,2	0,3	9

Name of the program	Description of the fields of application	t/min	Keyword Indication		
			1*	2*	3*
Relaxation Programs					
Rest 1	Program for deep, inner relaxation	30		X	X
Rest 2	For any present time stress	30		X	X
Sleep 1	Supportive in problems of falling or staying asleep	30		X	
Sleep 2	Regulates particularly the sleep rhythm	45		X	
Dream	Stimulating dream activity	45		X	
Conflict Resolution Programs					
Conflict Resolution Programs	Resolution and organization of individual conflicts	40	X		X
Children 1	Stress between ages 6 and 9	35	X		X
Children 2	Stress between ages 9 and 12	35	X		X
Gamma long	Resolution of blockages through the use of the meditative vibration	42			X
Gamma short	Resolution of blockages through specific reflex zones	10			X
Stress Programs					
Stress Basic	Resolution of psychological tension	30			X
Stress Immune	Regulation and enhancement of the immune system	35			X
Stress Hormone	Supports the regulation of hormonal symptoms	40			X
Stress Spasm	Basic program for support in cases of migraine and headaches	42			X
Depression/Psyche					
Psyche 1	Fatigue, weariness, listlessness	30	X		
Psyche 2	In the manic phase of depression	30	X		
Psyche 3	For endocrine depression	38	X		
Cerebral/Mental Programs					
Cerebral	Stimulation of brain activity, increasing wakefulness	35	X		X
Learning	Supports in case of learning disabilities and concentration problems	35			X
Memory	Increase of intellectual ability and stimulation of creativity	35			X
Additional Programs					
Waking	Exhaustion, reconvalescence, chronic fatigue, complaints of the elderly	35	X	X	
Degeneration	Basic therapy in all degenerative diseases	45	X		
Power-Nap	Regeneration program for the elimination of daily stress	15			X
Addiction	Basic program for all types of addiction	45	X		
*) 1 = Depression, 2 = Insomnia, 3 = Mental exhaustion/appearances of stress					

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esogetics GmbH
esogetics GmbH – Niederlassung